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Mix 12 grams (0.05 mole) of cobalt chloride hexahydrate, CoCl<sub>2</sub>·6H<sub>2</sub>O (237.93), and 8 grams (0.15 mole) of ammonium chloride (NH<sub>4</sub>Cl) with 10~20ml of water, and stir until most of the salts have dissolved; then add 25mL of concentrated ammonia solution and 0.2~0.3 gram of activated carbon, preferably a good grade taken from a freshly opened bottle. \*[Put the solution into a 100-mL Büchner flask which is fitted with a one-hole rubber stopper carrying a wide glass (10mm bore) reaching down to the bottom of flask. With the aspirator pump attached to the side arm of the flask, draw air through the solution at a brisk rate (Be careful that the water don't suck back from the aspirator in the event that the air inlet becomes clogged with crystals).]

After an hour or an hour and a half, the color of the solution in the flask should have changed from a reddish color to a yellowish brown. The color may easily be observed by swirling the liquid round the sides of the flask while holding against the light. Pass air for 10 minutes after all the red color has gone from the liquid; then filter the solution on a Büchner funnel. The product, mixed with carbon, remains on the filter paper.

Stir the filter cake into 100mL of water to which has been added 1mL of concentrated hydrochloric acid. This should be enough to give the solution a slight acid reaction after all the solid has been added. Test the solution with litmus; if necessary, add a few more drops of hydrochloric acid. Then heat to 50-60°C and filter while hot. The desired product is in the solution. It may be precipitated by using the common-ion effect. Add 15mL of concentrated hydrochloric acid to the hot solution and set aside to cool slowly. When it has cooled to room temperature, set the beaker in cool water or crushed ice. Filter the fine yellow crystals on a Büchner funnel; then wash, first with 30mL of 60 per cent alcohol and then with 30mL of 95 per cent alcohol. Suck as dry as possible; then dry on a watch glass in the oven at 80°C. The yield should be about 85 per cent.

The carbon catalyst used in this preparation not only catalyzes the oxidation but also catalyzes the attainment of equilibrium

$$[Co(NH_3)_5Cl]^{2+} + NH_3 \rightarrow [Co(NH_3)_6]^{3+} + Cl^{-}$$

which under the conditions of the reaction is greatly in favor of the desired [Co(NH<sub>3</sub>)<sub>6</sub>]<sup>3+</sup>.

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